

BEAUTIFUL SMILES BEGIN HERE

ENGINEERED AND OPTIMIZED
FOR AESTHETIC DENTAL
PROCESSES



KEROX Farben KEROX Effekten Marker

The perfect products for the combination of dipping and brushing technology for the special characterization of dental zirconia restorations.

Adapted for
KEROX HS, HT and UHT
zirconia materials

PREMIUM EUROPEAN QUALITY

- KEROX Farben
Colouring Liquid
- KEROX Effekten
- Marker
- Modelling Liquid

KEROX Farben Colouring Liquid



- 16 shades for the entire chromatic spectrum of the VITA classical A1-D6
- **Optimized for KEROX (HS, HT, UHT) materials**
- Water based with trace amounts of acid
- alcohol free so it doesn't evaporate
- Absorbs deep and doesn't leave spots on the surface
- Optimal dip time is 30 seconds
- Can be used in any sintering furnace

KEROX Effekten



- Use before sintering
- Available colours:
 - **Light Pink, Dark Pink**
 - Cervical & Fissures: **Brown**
 - Translucency: **Grey**

Marker



- Multi-brushing and multilayering
- Obtainable colours:
 - **Red**
 - **Blue**
 - **Brown**
 - **Green**

Modelling Liquid



KEROX Modelling Fluid is used to wet the dental technician's instruments during modelling.

It is also the perfect modelling aid for wetting the brush to disperse the material.
200ml

