

# ZIRCOSTAR MULTILAYER

HTML & UTML



HTML HIGH TRANSLUCENT MULTILAYER	UTML ULTRA HIGH TRANSLUCENT MULTILAYER
UNDEFINED LAYERS – MORE HOMOGENEOUS GRADIENTS	
COMPATIBLE WITH 71, 95 AND 98 SYSTEMS	
AVAILABLE IN 14MM AND 20MM THICKNESS	
UP TO 1389 MPa	UP TO 746 MPa
<b>APPLICATIONS:</b> FULL CONTOUR CROWNS & BRIDGES UP TO 3 UNITS	<b>APPLICATIONS:</b> FULL CONTOUR CROWNS & BRIDGES UP TO 3 UNITS FOR ANTERIOR RESTORATIONS

**PREMIUM  
EUROPEAN  
QUALITY**

AVAILABLE SHADES



A1    A2    A3    B2    D2

# KEROX DENTAL

35 YEARS OF EXPERIENCE  
IN HIGH PRECISION  
CERAMICS



## WHY KEROX MULTILAYER?

- High translucency, for life-like aesthetics
- High strength, for unparalleled mechanical properties
- Undefined layers → more homogeneous gradients
- Perfect Vita® Match



World Class Customer Support

## KEROX Ltd.

Kerox St. 1, Sósút, 2038 Hungary-EU

Phone: +36 23 560 700 EXT155

Fax: +36 23 545 215

E-mail: [info@keroxdental.net](mailto:info@keroxdental.net)