

HIGH TRANSLUCENT (HT)

FOR FULL CONTOUR &
FRAMEWORKS



OPTIMAL MILLING PROPERTIES FOR ALL
KIND OF RESTORATIONS

WHITE AND 5 SHADES



COMPATIBLE WITH ALL SYSTEMS
(71, 95, 98, 100 MM)
FLEXURAL STRENGTH UP TO 1351 MPa

APPLICATIONS:
FULL CONTOUR, FRAMEWORKS, INLAYS,
ONLAYS, CROWNS, COPINGS

**35 YEARS OF EXPERIENCE
IN HIGH PRECISION CERAMICS**

PREMIUM EUROPEAN QUALITY

WHY KEROX HT?

- Ideal for cut-back technology
- One of the highest fracture toughnesses on the dental market
- Perfect combination of density, strength and millability
- Proven longer lifetime for your milling burrs



IDEAL FOR CUT-BACK TECHNOLOGY

- Globally unique pressing technology → evenly more dense material composition → softer and lighter milling → less chipping
- The most widely practicable zirconia type for dental zirconia restorations
- More realistic and lifelike layering with porcelain for highly aesthetic restorations
- Perfectly smooth milling surface → less adjustments needed in the lab



World Class Customer Support

KEROX Ltd.

Kerox St. 1, Sósút, 2038 Hungary-EU

Phone: +36 23 560 700 EXT155

Fax: +36 23 545 215

E-mail: info@keroxdental.net