

KEROX PREMIA PMMA REFERENCES

Example: PA39818M

Kerox Premia	VITA	Blank	Blank	Blank
PMMA Blank	Shade	Diameter	Thickness	ML ?
1	2	3	4	5
P	A3	98	18	M

1	P	for all Kerox Premia PMMA blanks			
2 color of the blank	A1	A1 Vita shade			
	A2	A2 Vita shade			
	A3	A3 Vita shade			
	B1	B1 Vita shade			
	B2	B2 Vita shade			
	B3	B3 Vita shade			
	CL	Clear			
	CLO	Orange			
	CLB	Blue			
	CLBE	Berry			
	CLG	Green			
	PI	PINK PMMA			
	PIA1	Pink & A1 Vita shade 2 layers			
	PIA2	Pink & A2 Vita shade 2 layers			
	PIA3	Pink & A3 Vita shade 2 layers			
	BL1	Bleach 1 shade			
	BL2	Bleach 2 shade			
BL3	Bleach 3 shade				
4 Blank diameter (mm)	95	Diameter of the blank (mm)			
	71	Diameter of the blank (mm)			
	98	Diameter of the blank (mm)			
	115	Diameter of the blank (mm)			
5 Thickness (mm)	12, 14, 15, 16, 18, 20, 25,				
	30	Thickness of the blank (mm)			
6 (ML?)	M	Additional information if Multilayer PMMA			