

# ZIRCOSTAR REFERENCIA

Example: C00N9812AP

ZIRCOSTAR	VITA	Blank	Blank	Blank	Blank
Blank	Shade	material	diameter	thickness	Form
1	2	3	4	5	6
C	00	N	98	12	AP

1	C	for all ZIRCOSTAR blanks			
2 shade	00	White			
	AL	For light A-B Vita shades (HS)			
	AD	For dark A-B Vita shades (HS)			
	CL	For light C-D Vita shades (HS)			
	CD	For dark C-D Vita shades (HS)			
	A1	VITA shade (HT, UHT, HTML,UTML, 3D ML)			
	A2	VITA shade (HT, UHT, HTML,UTML, 3D ML)			
	A3	VITA shade (HT, UHT, HTML,UTML, 3D ML)			
	B2	VITA shade (HT, UHT, HTML,UTML, 3D ML)			
	D2	VITA shade (HT, UHT, HTML,UTML, 3D ML)			
	U1	Bleach 1 shade (HT)			
	U3	Bleach 1 shade (UHT)			
	U5	B1 shade (3D ML)			
	L1	Bleach 1 shade (3D ML)			
L2	Bleach 2 shade (3D ML)				
3 Blank material	N	HS	High strength		
	C	HTML	High Translucent		
	U	UHT	Ultra High Translucent		
	V	HTML	High Translucent Multilayer		
	W	UTML	Ultra High Translucent Multilayer		
	Z	3D ML	3D Multilayer		
4 Blank diameter	95	diameter of the blank (mm)			
	71	diameter of the blank (mm)			
	98	diameter of the blank (mm)			
	100	diameter of the blank (mm)			
5 Blank thickness	10,12,14, 16,18,20,2 5	mm thick blank			
6 Blank form	ZZ	Zirkon Zahn compatible			
	AG	Amann Girrbach compatible			
	AP	Wieland compatible			
	DM	Dentmill compatible			